



Numedico's Medical Device Product Range 2023

The Numedico medical device product range is a collection of standard and safety syringe products which are ideal for hospitals, healthcare clinics and at home health care.



Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com



About Numedico

Focused on safety.
Driven by excellence.
Doing things differently.

Established in 2016 out of a desire to pursue optimal medtech solutions to challenge the status quo, Numedico is a globally focused medical supplies company based in Australia.

Our solutions include safety syringes for the use at home, hospitals and clinics.

We care about our customers and with offices in various global hubs are able provide an excellent customer service and a hassle-free experience.



Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com

Table of Contents

Numedico's Medical Device Product Range



**ClickZip™ Needle Retractable
Safety Syringe Range**

4



Safety Syringe Range

7



Standard Syringe Range

9



Insulin Syringe Range

11



Safety Blood Collection Range

13

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

PREMIUM RANGE

ClickZip™ Manual Retractable



Description of Technology & Product Features



ClickZip™ prevents reuse and reduces needlestick injury due to its innovative tilted needle technology that encapsulates the needle in the syringe barrel.

The innovative ClickZip™ range includes the Manual Retractable Safety Insulin Syringe with needle and the ClickZip™ Needle Retractable Safety Syringe featuring a high-quality, triple bevel needle designed for increased patient comfort. The unique design also ensures minimal leakage and medicine wastage.

Ideal for:



Hospitals



Home Health Care



Healthcare Clinic



Veterinary



Cosmetic

Features:

Material

- Highest quality Japanese stainless steel cannulas
- Latex free

Technology

- Retraction mechanism, which tilts the needle once retracted
- Luer lock needle makes it exchangeable to draw medication
- Reuse prevention (RUP)
- Minimal dead space of <.02ml

Functionality

- No change to the injection technique
- The user can hear and feel a click when the safety mechanism is activated
- The user manually retracts the needle with the plunger
- After retraction the needle tilts to one side in the barrel preventing reuse
- Meets NIOSH recommendation for needlestick devices.



All products have 510K FDA approval.

Providing quality medical supplies globally from a trusted Australian brand.


www.numedico.com

www.clickzipsafety.com



Product Structure of the ClickZip™ technology

ClickZip™ Needle Retractable Safety Syringe is a single use sterile needle retractable safety syringe (Reuse Prevention Syringe) that is designed to draw the needle up into the syringe barrel to reduce needlestick injury and prevent reuse.



ClickZip™ is the future of needle safety.

All our syringes meet the necessary ISO standards.



Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com

Order Information for ClickZip Range



Description	Needle Gauge			Needle Length
40U ClickZip Manual Retractable Safety Insulin Syringe with needle	29	30		0.5"
100U ClickZip Manual Retractable Safety Insulin Syringe with needle	29	30		0.5"
0.5ml ClickZip Manual Retractable Safety Syringe with needle	23	25	28	5/8" 1"
1ml ClickZip Manual Retractable Safety Syringe with needle	23	24	25	5/8" 1"
3ml ClickZip Needle Retractable Safety Syringe with needle	21	22	23 24	1" 1.5"
5ml ClickZip Needle Retractable Safety Syringe with needle	23	25		5/8" 1"
10ml ClickZip Needle Retractable Safety Syringe with needle	21	22	23 24	1"
20ml ClickZip Needle Retractable Safety Syringe with needle	23	24	25	1"
50ml ClickZip Needle Retractable Safety Syringe with needle	23	24	25	1"

For details on quantity per carton, order codes, or other SKU availability email: sales.usa@numedico.com



All products are stocked in the Numedico warehouses and ready for immediate delivery and have 510K FDA approval.
For order enquiries email: sales.usa@numedico.com

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com



STANDARD RANGE

Safety Cap Needle



Description of Technology & Product Features

Made using high-quality cannulas for increased patient comfort, they are available in a range of sizes and styles depending on clients needs, and can be delivered worldwide.

Ideal for:



Hospitals



Home Health Care



Healthcare Clinic



Veterinary



Cosmetic

Features:

Material:

- Made of Polypropylene/Polyisoprene
- Latex-free

Functionality:

- Safety cap to prevent needle stick injury
- Safe plunger backstop



All products have 510K FDA approval.

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com

Order Information for Safety Cap Needle Range

The products in this range include:

- Standard safety needle – luer lock | luer slip | universal
- Syringe with safety needle – luer lock | luer slip
- Syringe only – luer lock | luer slip

Description	Needle Gauge	Needle Length
Syringe Only 3ml Luer Lock No Needle	n/a	n/a
Syringe Only 3ml Luer Slip No Needle	n/a	n/a
Safety Needle Universal	23	1.5" 1"
Disposable Syringe 1ml with Needle Luer Lock	23	1"

For details on quantity per carton, order codes, or other SKU availability email: sales.usa@numedico.com



All products are stocked in the Numedico warehouses and ready for immediate delivery and have 510K FDA approval. For order enquiries email: sales.usa@numedico.com

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com

Standard Syringe Range



Description of Technology & Product Features

The Numedico Everyday Standard Syringe Range is comprised of a range of traditional standard and safety syringe products which are ideal for:



Hospitals



Home Health Care



Healthcare Clinic



Veterinary




Cosmetic

Made using high-quality cannulas for increased patient comfort, they are available in a range of sizes and styles depending on clients needs, and can be delivered worldwide. Features include:

- Hypodermic Syringes for single use with or without needle
- Made of Polypropylene/Polyisoprene
- Highly transparent barrel
- High contrast black graduation
- Safe plunger backstop
- Latex-free
- CE Marked
- US 510k Cleared
- ISO 13485:2003 Certified

Syringe Only



Syringe with Needle



Needle Only



All our syringes meet the necessary ISO standards.

Order Information for Standard Syringe Range

The products in this range include:

- Standard syringe with needle – luer lock | luer slip | fixed
- Standard needle – luer lock | luer slip
- Syringe only – luer lock | luer slip

Product	Sizes Available
Syringe Only	0.5ml-60ml
Needle Only	19G-27G 25mm-38mm
Syringe with Needle	Any combination from above

For more details and for full specifications contact: sales.usa@numedico.com

Insulin Syringe Range



Description of Technology & Product Features

Our diabetes syringe range offers our consumers a high quality insulin syringe for the comfortable administration of insulin.



Hospitals



Home Health Care



Healthcare Clinic

Product Features

- Available in various sizes to suit personal needs
- Clear barrel to ensure rapid, precise volume reading
- Ultra-fine needle to enhance patient comfort and reduce pain
- Latex free
- ClickZip™ - Needle retracts into the syringe barrel as a safety feature to prevent needlestick injuries
- Low dead space - improved accuracy and less medication waste when compared to conventional syringes
- Bold and clear graduated insulin unit markings for accurate dosing

Insulin Syringe Range



All our syringes meet the necessary ISO standards.

The products in this range include:

- Insulin syringe with fixed needle

Syringe Capacity	Needle		Packaging
	Gauge	Length	
0.3cc/mL 0.5cc/mL 1.0cc/mL	27G	1/2" (13mm) 5/16" (8mm)	100 pcs/inner box 2000 pcs/carton
	28G	1/2" (13mm) 5/16" (8mm)	
	29G	1/2" (13mm) 5/16" (8mm)	
	30G	1/2" (13mm) 5/16" (8mm)	
	31G	1/2" (13mm) 5/16" (8mm)	

For more details and for full specifications contact: sales@numedico.com.au

OTHER RANGE

Safety Blood Collection Sets



Description of Technology

Numedico Safety Blood Collection Needle is a safety-engineered device intended for venous blood collection. The safety mechanism can be activated immediately after blood draw, helping to protect against needlestick injury.

Ideal for:



Hospitals



Healthcare Clinic



Cosmetics



Dental

Product Features

- **Easy handling**
 - Makes venipuncture procedure more comfortable for patients and healthcare workers
 - Provides other indications for short IV infusion
 - Convenient for poor veins



All products have 510K FDA approval.

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com

Order Information for Safety Blood Collection Sets

Safety Blood Collection Needle without pre-attached holder

Description	Tubing Length	Needle Gauge		Needle Length
Safety Blood Collection Sets Without Holder	12"	23	25	3/4"

For details on quantity per carton, order codes, or other SKU availability email: sales.usa@numedico.com



All products are stocked in the Numedico warehouses and ready for immediate delivery and have 510K FDA approval. For order enquiries email: sales.usa@numedico.com

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com

www.clickzipsafety.com



Our future focus is driven by a desire to improve healthcare and to find the best and safest solutions to suit our customer's needs now and tomorrow.

For more details and for full specifications contact: **sales@numedico.com.au**

Providing quality medical supplies globally from a trusted Australian brand.

www.numedico.com.au