



Numedico
TECHNOLOGIES

Numedico Everyday Standard Syringe Range

The Numedico Everyday Standard Syringe Range is a range of traditional standard and safety syringe products which are ideal for hospitals, healthcare clinics, or for at home health care.



Numedico – providing quality medical supplies globally from a trusted Australian brand



Numedico
TECHNOLOGIES

Numedico Technologies

Focused on safety | Driven by excellence | Doing things differently

Established in 2016 out of a desire to pursue optimal medtech solutions and to challenge the status quo, Numedico is a globally focused safety medical devices company based in Australia.

Our future focus is driven by a desire to improve healthcare and to find the best and safest solutions to suit our customer's needs now and tomorrow. This includes safety in medical technology, in home health care innovations and alternative drug delivery options.

We push ourselves to be better today than we were yesterday, building business through best practice and collaborative partnerships. You can trust our team to be responsive, transparent, accountable, and honest. We care about our customers and with offices in various global hubs are able provide an excellent customer service and a hassle-free experience.

Our new look

In this new era of medical technologies, Numedico's reputation as a leader in the medical devices industry is strengthened by our existing and growing brand portfolio of innovative medical technology.


Description of technology

The Numedico Everyday Standard Syringe Range is comprised of a range of traditional standard and safety syringe products which are ideal for hospitals, healthcare clinics, home health care, veterinary and cosmetic.

Made using high-quality cannulas for increased patient comfort, they are available in a range of sizes and styles depending on clients needs, and can be delivered worldwide.

- Hypodermic Syringes for single use with or without needle
- Made of Polypropylene/Polysoprene
- Highly transparent barrel
- High contrast black graduation
- Safe plunger backstop
- Latex-free
- CE Marked
- US 510k Cleared
- ISO 13485:2003 Certified






We strive to consistently deliver the best-quality products to our client community.








Note: Packaging image is for illustrative purposes only, not actual package.

A view of the Numedico Everyday Standard and Safety Syringe Range

The Numedico Everyday standard and safety syringe range as shown below is available in the sizes and types shown in the table on the following page.

				
1ml Safety Needle Luer Lock	1ml Standard Syringe Luer Lock	3ml Luer Lock	10ml Luer slip	1ml Standard Syringe Luer Lock

				
1ml Standard Syringe Luer Slip	1ml Standard Syringe Luer Lock Side Needle	3ml Safety Needle Luer Lock	Safety Needle Luer Slip Lock	Blood collection needle

List of standards for Numedico Everyday Safety Syringe Range


The Numedico Everyday technology pertains to the international standard ISO7864, ISO7886, as the ISO standards that pertains to our range.

Numedico Everyday Safety and Standard Syringe Range

The products in this range include:

- Standard safety needle – luer lock | luer slip
- Standard needle – luer lock | luer slip | fixed
- Syringe only – luer lock | luer slip

Needle	Size mm	Size inch	Hub colour	Insulin 1ml	0.5ml	1ml	3ml	5ml	10ml
30G	13	1/2		✓	✓	✓			
29G	13	1/2		✓	✓	✓			
27G	8	5/16		✓	✓	✓			
26G	13	1/2		✓		✓			
26G	16	5/8		✓		✓			
26G	7.9	5/16		✓		✓			
25G	25	1				✓	✓	✓	✓
25G	16	5/8					✓	✓	✓
24G	25	1				✓	✓	✓	✓
24G	38	1 1/2					✓	✓	✓
24G	16	5/8					✓	✓	
23G	25	1				✓	✓	✓	✓
23G	38	1 1/2					✓	✓	✓
23G	32	1 1/4					✓	✓	✓
22G	25	1					✓	✓	✓
22G	38	1 1/2					✓	✓	✓
22G	32	1 1/4						✓	✓
21G	38	1 1/2							✓




Packaging inner box and carton quantities by syringe size

Needle Size	0.5ml	1ml	3ml	5ml	10ml
Qty Box/pcs	100	100	100	100	100
Qty Carton/pcs	3,000	3,000	2,400	1,800	1,200
Qty per container	2,000,000	1,560,000	1,104,000	1,000,000	1,000,000

Note:

- The product offers low dead space (exceeding standard requirement).
- Prevent heavy pressing, keep out direct sunlight, extreme heat, exposure to rain, snow or wetting.
- Keep at a temperature between +5°C and +60°C.
- Sterilization: EO gas.
- Shelf Life: 5 years from manufacturing date.
- Volume size may vary from the table above.



Numedico Technologies Pty Ltd
Registered office
Level 21, 91 King William Street SA 5000 AUSTRALIA
sales@numedico.com.au



Numedico
TECHNOLOGIES
numedico.com

Numedico – providing quality medical supplies globally from a trusted Australian brand